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Aim

Comparison of different antigen rapid tests using uniform sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs which are used as test specimens included in an evaluation panel.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable Cq values diluted 1:10 in negative samples in PBS (Cq = Quantification cycle; also reported as CT = cycle threshold). The Cq values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standard. In the case of the PCRs used, a Cq value of 25 corresponds to around 10⁶ RNA copies/mL. The replication of the virus in cell culture as another characteristic of the samples was determined as a possible correlate for infectiousness.

Due to the high number of tests under evaluation few specimens of the first panel were replaced (resulting in panel 1 versions 1 and 2), and a new panel (panel 2) was manufactured in analogous manner. For panel 1 versions 1 and 2, 18 or 17 samples, respectively, were characterized as Cq ≤25, 23 samples as Cq 25-30, and 9 or 10 samples, respectively, as Cq ≥30; the composition of panel 2 is: 20 samples with Cq ≤25, 20 samples with Cq 25-30, and 10 samples with Cq ≥30.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50µl of the pool (panel 1 versions 1 and 2), or 25µl (panel 2) were analysed using the components of the test provided, e.g. swabs and buffer. In case of tests designed for various sample types (e.g. nasopharyngeal, nasal, saliva, etc), the test specific swab dedicated for

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Kontakt:

E-Mail: sarscov2ivd@pei.de

nasopharyngeal samples was used. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers using the uniform sample sets allows an overview of the current state of the art regarding sensitivity.

The current state of the art was defined as corresponding to a minimal sensitivity of 75% for the pools with $Cq \leq 25$. The results do not allow any conclusions regarding specificity of the tests.

Table 1 summarizes the results of those POCTs assessed as reflecting the current state of the art (alphabetical order of manufacturers).

Table 2 summarizes the results of those POCT (alphabetical order of manufacturers) which do not fulfill the sensitivity criteria and consequently are deleted from the list of reimbursable tests (see the homepage of Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM)).

The first part of this evaluation has already been published (1); this publication summarizes further details of the comparative evaluation.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for reimbursement, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

Contact

Email: sarscov2ivd@pei

(1) Scheiblauer Heinrich, Filomena Angela, Nitsche Andreas, Puyskens Andreas, Corman Victor M, Drosten Christian, Zwirgmaier Karin, Lange Constanze, Emmerich Petra, Müller Michael, Knauer Olivia, Nübling C Micha. Comparative sensitivity evaluation for 122 CE-marked rapid diagnostic tests for SARS-CoV-2 antigen, Germany, September 2020 to April 2021. Euro Surveill. 2021;26(44):pii=2100441. <https://doi.org/10.2807/1560-7917.ES.2021.26.44.2100441>

Sensitivität / Sensitivity

| AT-Nr. / AT-No. | AT-Nr. / Selbsttest / AT-No. self-test | Ref-Nr. / ID-No. * | Hersteller / Manufacturer | Testname / Test name | Zielantigen / target antigen | Cq ≤25 | Cq 25-30 | Cq ≥30 | Gesamt-Sensitivität / total sensitivity |
|-----------------|--|--------------------|---|--|------------------------------|--------|----------|--------|---|
| AT022/20 | | F-NCOV-01G | SD Biosensor, Inc. | STANDARD F COVID-19 Ag FIA | N | 100,0% | 65,2% | 0,0% | 66,0% |
| AT1064/21 | | 2.13.1.587 | Senova Gesellschaft für Biowissenschaft und Technik mbH | GreenLight SARS-CoV-2 Antigen-Test | N | 95,0% | 10,0% | 0,0% | 42,0% |
| AT918/21 | AT1291/21 / AT1338/21 | SCVC02 | SGA Mühendislik DAN. EG. love DIS.Ltd.STI | V-Chek SARS-CoV-2 Rapid Ag Test Kit (Colloidal Gold) | N | 94,1% | 26,1% | 0,0% | 44,0% |
| AT748/21 | | n.a. | Shanghai Liangrun Biomedicine Technology Co., Ltd. | LionRun™ SARS-CoV-2 Antigen Rapid Test Kit (Lateral Flow Immunoassay) | N | 100,0% | 85,0% | 0,0% | 74,0% |
| AT857/21 | | n.a. | Shen Zhen Zi Jian Biotechnology Co., Ltd. | Zijian SARS-CoV-2 (COVID-19) Antigen Rapid Detection Kit (Lateral Flow Method) | N | 80,0% | 10,0% | 0,0% | 36,0% |
| AT928/21 | | CG201A | Shenzhen CAS-Envision Medical Technology Co., Ltd. | SARS-CoV-2 Antigen Rapid Detection Kit | N | 100,0% | 100,0% | 70,0% | 94,0% |
| AT260/21 | | n.a. | Shenzhen Dymind Biotechnology Co., Ltd. | SARS-CoV-2 Antigen Test Kit (Colloidal Gold) | N | 100,0% | 30,0% | 0,0% | 52,0% |
| AT829/21 | | n.a. | Shenzhen Dymind Biotechnology Co., Ltd. | SARS-CoV-2 Antigen Test Kit (Colloidal Gold) | N | 100,0% | 30,0% | 0,0% | 52,0% |
| AT458/21 | | HRK-6820 | Shenzhen Huaree Technology Co., Ltd. | SARS-CoV-2 Antigen Saliva Rapid Test Kit | N | 100,0% | 90,0% | 10,0% | 78,0% |
| AT424/21 | | HRK-6620 | Shenzhen Huaree Technology Co., Ltd. | SARS-CoV-2 Antigen Rapid Test Kit (Immunochromatography) | N | 100,0% | 45,0% | 0,0% | 58,0% |
| AT472/21 | | 203-020 | Shenzhen Huian Biosci Technology Co., Ltd. | SARS-CoV-2 Antigen Test Kit (Colloidal Gold) | N | 100,0% | 60,0% | 0,0% | 64,0% |
| AT1065/21 | | 808151-1 | Shenzhen Kang Sheng Bao Bio- Technology Co., Ltd. | SARS-CoV-2 Antigen Rapid Test (Lateral Flow Assay) | N | 95,0% | 10,0% | 0,0% | 42,0% |
| AT1019/21 | | 150129 | Shenzhen Kingfocus Biomedical Engineering Co., Ltd. | COVID-19 Antigen Detection Kit (Quantum Dots-Based Immunofluorescence Chromatography) | N | 100,0% | 95,0% | 60,0% | 90,0% |
| AT420/21 | | K602-20 | Shenzhen Kisshealth Biotechnology Co., Ltd. | KISSH SARS-CoV-2 Antigen Test Kit (GICA) | N | 100,0% | 85,0% | 10,0% | 76,0% |
| AT417/20 | | n.a. | Shenzhen Lvshiyuan Biotechnology Co., Ltd. | Green Spring SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) | N | 100,0% | 95,7% | 40,0% | 86,0% |
| AT070/20 | | MF-60 | Shenzhen Microprofit Biotech Co., Ltd. | fluorecare COVID-19 SARS-CoV-2 Spike Protein Test Kit (Colloidal Gold Chromatographic Immunoassay) | S | 100,0% | 47,8% | 10,0% | 58,0% |
| AT071/20 | | MF-63 | Shenzhen Microprofit Biotech Co., Ltd. | SARS-CoV-2 Spike Protein Test Kit (Fluorescence Immunoassay) | S | 100,0% | 85,0% | 10,0% | 76,0% |
| AT544/21 | AT1308/21 | MF-68 | Shenzhen Microprofit Biotech Co., Ltd. | fluorecare SARS-CoV-2 Antigen Test Kit | N | 90,0% | 15,0% | 0,0% | 42,0% |
| AT844/21 | | RNS92048B | Shenzhen Reagent Technology Co., Ltd. | SARS-CoV-2 Antigen IVD kit SWAB | N+S | 100,0% | 50,0% | 0,0% | 60,0% |
| AT108/20 | | LFA0401+25N | Shenzhen Watmind Medical Co., Ltd. | SARS-CoV-2 Ag Diagnostic Test Kit (Colloidal Gold) | N | 100,0% | 95,7% | 20,0% | 82,0% |
| AT385/20 | | LFB0401-25N | Shenzhen Watmind Medical Co., Ltd. | SARS-CoV-2 Ag Diagnostic Test Kit (Immunofluorescence) | N | 100,0% | 60,9% | 0,0% | 62,0% |
| AT409/21 | | G86247 | Shenzhen YHLO Biotech Co., Ltd. | GLINE-2019-nCoV Ag | N | 90,0% | 30,0% | 0,0% | 48,0% |